

Technical Information

Adiponectin, Mouse

Total Adiponectin

Cat. No.:	E091M
Tests:	96
Method:	ELISA
Range:	0.025 - 1 ng/ml native Adiponectin
Sensitivity:	~ 0.01 ng/ml
Incubation time:	3 hours
Sample volume:	100 µl (after dilution 1:10'000)
Sample type:	Serum and plasma
Sample preparation:	Generally, samples should be refrigerated as soon as possible following collection. Samples are stable maximum 2 days at room temperature. Long-term storage up to 2 years at -20 °C or below. Maximum 5 freeze and thaw cycles.

Species: Mouse

Intended use:

Adiponectin is a 30kDa protein and its percentage in serum proteins is 0.01 %. In vivo, it appears with different oligomers and it is mainly synthesized by adipocytes. Until now, IGF-1 is the only known natural inductor of synthesis.

Low Adiponectin levels are closely associated with insulin resistance and metabolic syndrome as well as an increased risk of type 2 diabetes mellitus and cardiovascular disease.

Today, Adiponectin is thought to act as an endogenous insulin sensitizer by decreasing excessive glucose levels without increasing insulin concentrations and by stimulating the burning of adipose tissue in muscle and liver.

Adiponectin is associated with glucose and lipid metabolism and is assumed to have direct antiatherogenic characteristics. Furthermore, it is involved in inflammatory processes.

Clinical significance:

- Obesity
- Arteriosclerosis
- Energy metabolism
- Coronary diseases

For further information please contact / Für weitere Informationen wenden Sie sich bitte an /
Pour plus d'informations, veuillez contacter:

www.tecomedical.com



A EUROBIO SCIENTIFIC COMPANY

Switzerland / Headquarters

TECOmedical AG
Gewerbstrasse 10
4450 Sissach
Phone +41 61 985 81 00
Fax +41 61 985 81 09
Mail info@tecomedical.com

Germany

TECOmedical GmbH
Wasserbreite 57
32257 Bünde
Phone +49 52 23 985 99 99
Fax +49 52 23 985 99 98
Mail info@tecomedical.com

Benelux

TECOmedical Benelux BV
Prins Willem-Alexanderlaan 301
7311 SW Apeldoorn, The Netherlands
Phone +31 30 307 87 30
Fax +31 30 307 49 39
Mail benelux@tecomedical.com

Austria

TECOmedical AG
Phone 0800 20 40 66
Fax 0800 20 40 55
Mail info@tecomedical.com