

## Information

*always your partner*

### TECO<sup>®</sup> Fast Candida Ag & IgG



2 lateral flow assays for rapid, parallel diagnosis of invasive candidiasis and candidemia on the antigen and antibody level

- Fluorescence-based lateral flow assays designed for readout via reader
- (Semi-)quantitative, calculation of IgG concentration in AU/mL and Mannan Antigen in pg/mL
- Validated for serum samples
- 2 results in 40 minutes obtained from 305  $\mu$ L of serum
- Detection limit IgG: 7 AU/mL, cut-off value of 135 AU/mL
- Detection limit Antigen: 21.5 pg/mL, cut-off value of 60 pg/mL
- Increase in diagnostic significance through combination of two parameters

	<b>Fast Candida Mannan IgG</b>	<b>Fast Candida Mannan Antigen</b>
<b>Clinical sensitivity</b>	86.2% (n = 222)	88.1% (n = 155)
<b>Clinical specificity</b>	95.3% (n = 222)	96.9% (n = 155)



## TECO® *Fast* Candida IgG and Antigen LFA – sample preparation



TECO® *Fast* Candida Mannan Antigen LFA  
Catalogue No.: TE1081

&

TECO® *Fast* Candida IgG Antibody LFA  
Catalogue No.: TE1083

Für weitere Informationen wenden Sie sich bitte an:  
For further information please contact:  
Pour plus d'information, veuillez contacter:

[www.tecomedical.com](http://www.tecomedical.com)



A EUROBIO SCIENTIFIC COMPANY

Switzerland / Headquarters  
**TECOmedical AG**  
Gewerbstrasse 10  
4450 Sissach  
Phone +41 61 985 81 00  
Fax +41 61 985 81 09  
Mail [info@tecomedical.com](mailto:info@tecomedical.com)

Germany  
**TECOmedical GmbH**  
Wasserbreite 57  
32257 Bünde  
Phone +49 52 23 985 99 99  
Fax +49 52 23 985 99 98  
Mail [info@tecomedical.com](mailto:info@tecomedical.com)

Austria  
**TECOmedical AG**  
Phone 0800 20 40 66  
Fax 0800 20 40 55  
Mail [info@tecomedical.com](mailto:info@tecomedical.com)